




*Precision
Crafted.
Excellence
Delivered.*



www.venhawk.com

SCALPELS



Precision Engineered. Trusted By Professionals.

Dissecting Knives

Anatomische Skalpelle
Bisturíes para disección
Bistouris à dissection
Bisturi per dissezione

Scalpels

Skalpelle
Bisturíes
Bistouris
Bisturi



Brain Knives

Hirnmesser

Cuchillos para cerebro

Couteaux à cerveau

Coltelli per il cervello



V2



1/2

Virchow

V2-003-16

V2-003-24

with hollow handle
mit Hohlheft
con mango hueco
avec manche creux
con manico vuoto




160 mm




200 mm




240 mm



V2-003-16
26 cm / 10 1/4"



V2-003-20
33.5 cm / 13 1/4"



V2-003-24
38 cm / 15"

Blade Holders / Blade Breakers

Klingenhalter / Klingenbrecher

Porta-hojas / Rompe-hojas de afeitar

Pincettes pour saisir et casser les lames de rasoir

Pinze per afferrare e rompere le lamette di rasoio


Scalpels

Skalpelle

Bisturíes

Bistouris

Bisturi



handle diameter
Griffdurchmesser
diámetro del mango
diamètre du manche
diametro di manico
5,5 mm

Troutman-Chris

V2-004-09

9 cm / 3 4/8"




handle diameter
Griffdurchmesser
diámetro del mango
diamètre du manche
diametro di manico
5,5 mm

Troutman-Chris

V2-005-09


9 cm / 3 4/8"



1/1

V2-005-11

11 cm / 4 3/8"



handle diameter
Griffdurchmesser
diámetro del mango
diamètre du manche
diametro di manico
8 mm

Castroviejo

V2-006-13

13 cm / 5 1/8"

Scalpel Handles

Skalpelli-griffe

Mangos de bisturí

Manches pour bistouris

Manici per bisturi



V2



1/2 Nr. 4

V2-007-04

13.5 cm / 5 3/8"

solid
massiv
macizo
solide
pieno



Fig. 18-36



1/2 Nr. 3 & 4

V2-007-34

16 cm / 6 2/8"

solid
massiv
macizo
solide
pieno



Fig. 10-36



1/2 Nr. 4L

V2-008-41

21 cm / 8 2/8"

solid
massiv
macizo
solide
pieno



Fig. 18-36



1/2 Nr. 4L

V2-008-45

20.5 cm / 8 1/8"

solid
massiv
macizo
solide
pieno



Fig. 18-36

Blade Holders / Blade Breakers

Klingenhalter / Klingenbrecher

Porta-hojas / Rompe-hojas de afeitár

Pincettes pour saisir et casser les lames de rasoir

Pinze per afferrare e rompere le lamette di rasoio

Scalpels

Skalpelle

Bisturís

Bistouris

Bisturi



Collin

V2-009-08

13 cm / 5 3/8"

with hollow handle
mit Hohlheft

con mango hueco
avec manche creux
con manico vuoto



Collin

V2-010-09

14.5 cm / 5 6/8"

with hollow handle
mit Hohlheft

con mango hueco
avec manche creux
con manico vuoto



V2-011-04

13.5 cm / 5 3/8"

with plastic handle
mit Kunststoffheft

con mango plástico
avec manche en matière plastique
con mango plástico



V2-012-03

15.5 cm / 6 1/8"

swivel head
Schwenkkopf
con cabeza orientable
à tête pivotante
con testa inclinabile



V2-013-04

16.5 cm / 6 2/8"

with hollow handle
mit Hohlheft
con mango hueco
avec manche creux
con manico vuoto

Scalpel Handles

Skalpellgriffe

Mangos de bisturí

Manches pour bistouris

Manici per bisturi



V2



Fig. 10-17

1/2 Nr.3

V2-014-03

12 cm / 4 6/8"

solid
massiv
macizo
solide
pieno



Fig. 10-17

1/2 Nr.3

V2-014-04

12 cm / 4 6/8"

solid
massiv
macizo
solide
pieno



Fig. 10-17

1/2 Nr.3

V2-014-05

12 cm / 4 6/8"

solid, cm-graduation
massiv, mit Zentimeter-
Graduierung
macizo, graduado en cm
solide, gradué en cm
pieno, graduato in cm



Fig. 10-17

1/2 Nr.3L

V2-015-31

21.5 cm / 8 4/8"

solid, long
massiv, lang
macizo, largo
solide, long
pieno, lungo




Fig. 10-17

1/2 Nr.3L

V2-015-35

20.5 cm / 8 1/8"

solid, angled, long
massiv, abgewinkelt, lang
macizo, angulado, largo
solide, coudé, long
pieno, angolare, lungo

Scalpel Handles

Skalpellgriffe
Mangos de bisturí
Manches pour bistouris
Manici per bisturi

Scalpels

Skalpelle
Bisturíes
Bistouris
Bisturi



Fig. 10-17

$\frac{1}{2}$ Nr.7

V2-016-07

16 cm / 6 $\frac{2}{8}$ "

solid
massiv
macizo
solide
pieno



Fig. 10-17

$\frac{1}{2}$ Nr.7K

V2-016-71

12.5 cm / 4 $\frac{7}{8}$ "

solid
massiv
macizo
solide
pieno



Fig. 10-17

$\frac{1}{2}$ Nr.3

V2-017-03

15.5 cm / 6 $\frac{1}{8}$ "

with round hollow handle
mit rundem Hohlheft
con mango hueco redondo
avec manche creux rond
con manico vuoto rotondo



Fig. 10-17

$\frac{1}{2}$ Nr.3

V2-018-03

12.5 cm / 4 $\frac{7}{8}$ "

with plastic handle
mit Kunststoffheft
con mango plástico
avec manche en matière plastique
con manico di materia plastica

Scalpel Handles

Skalpellgriffe

Mangos de bisturí

Manches pour bistouris

Manici per bisturi



V2



1/2 Nr.3



Fig. 10-17

V2-019-03

14 cm / 5 4/8"

solid
massiv
macizo
solide
pieno



1/2 Nr.3



Fig. 10-17

V2-020-03

14 cm / 5 4/8"

solid
massiv
macizo
solide
pieno



1/2 Nr.3



Fig. 10-17

Kaye V2-021-03

15 cm / 5 7/8"

solid
massiv
macizo
solide
pieno



1/2

V2-022-13

13.5 cm / 5 3/8"

for micro scalpel blades
für Mikro-Skalpellklingen
para micro hojas
pour micro-lames de bisturis
per micro lame

Amputation Knives

Amputationsmesser
Cuchillos para amputación
Couteaux à amputation
Coltelli per amputazione

Scalpels

Skalpelle
Bisturíes
Bistouris
Bisturi



1/2

Liston V2-023-16

29 cm / 11 ³/₈"
with hollow handle
mit Hohlheft
con mango hueco
avec manche creux
con manico vuoto



1/2

Liston V2-023-19

32 cm / 12 ⁵/₈"
with hollow handle
mit Hohlheft
con mango hueco
avec manche creux
con manico vuoto

Resection Knives

Resektionsmesser

Cuchillos para resección

Couteaux à résection

Coltelli per resezione



V2



V2-024-06

19 cm / 7 ⁴/₈ "

solid handle, with dissecting end
massives Heft mit Schaber-Ende
mango macizo, con periostótomo
manche solide à rugine
con manico pieno con raschino



Langenbeck V2-025-05

18 cm / 7 ¹/₈ "

with hollow handle
mit Hohlheft
con mango hueco
avec manche creux
con manico vuoto



Langenbeck V2-026-05

18 cm / 7 ¹/₈ "

with hollow handle
mit Hohlheft
con mango hueco
avec manche creux
con manico vuoto



Dermatome

Handdermatom

Dermatomo

Dermatome

Dermatomo a mano

Silver V2-027-19

19.5 cm / 7 ⁵/₈ "

for usual razor blades
für handelsübliche Rasierklingen
para hojas de afeitar usuales
pour lames de rasoir modèle standard
per lamette da barba normali

Dermatomes

Dermatome
Dermátomos
Dermatomes
Dermatomi

Scalpels

Skalpelle
Bisturíes
Bistouris
Bisturi



1/2

Humby V2-028-00

32 cm / 12 5/8"
complete with blade
komplett mit Klinge
completo con hoja
complet avec lame
completo con lama



1/2

Humby V2-028-01

Spare Blade
Ersatzklinge
Hoja de repuesto
Lame de rechange, seule
Lama di ricambio sola



V2-028-00



1/2

Schink V2-029-00

30 cm / 11 6/8"
adjustable from 0.1 to 2.0 mm
einstellbar von 0,1 bis 2,0 mm
ajustable de 0,1 hasta 2,0 mm
réglable en épaisseur de 0,1 à 2,0 mm
regolabile da 1,0 a 2,0 mm

complete with 2 skin fixing plates and one blade
komplett mit 2 Spannplatten und einer Klinge
complet avec deux plaques pour tendre la peau et une lame
completo con 2 placas para extender la piel y una hoja
completo con 2 placche per stendere la pelle ed una lama



1/2



V2-029-00

Schink V2-029-01

Ersatzklinge
Spare Blade
Hoja de repuesto
Lame de rechange, seule
Lama di ricambio sola



V2-029-00



1/2

Schink V2-029-02

Skin Fixing Plate
Spannplatte
Placa para extender la piel
Plaque pour tendre la peau
Placca per stendere la pelle